

Product datasheet for **GA102373**

Human IL1 Receptor I (IL1R1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL1R1 CRISPRa kit - CRISPR gene activation of human interleukin 1 receptor type 1
Symbol:	IL1R1
UniProt ID:	P14778
Synonyms:	CD121A; D2S1473; IL-1R-alpha; IL1R; IL1RA; P80
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)