

Product datasheet for **GA102374**

Human IL1RAP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL1RAP CRISPRa kit - CRISPR gene activation of human interleukin 1 receptor accessory protein
Symbol:	IL1RAP
UniProt ID:	<u>Q9NPH3</u>
Synonyms:	C3orf13; IL-1RAcP; IL1R3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)