

Product datasheet for **GA102381**

Human IL3RA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL3RA CRISPRa kit - CRISPR gene activation of human interleukin 3 receptor subunit alpha
Symbol:	IL3RA
UniProt ID:	P26951
Synonyms:	CD123; hIL-3Ra; IL3R; IL3RAY; IL3RX; IL3RY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)