

Product datasheet for **GA102383**

Human IL4R activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL4R CRISPRa kit - CRISPR gene activation of human interleukin 4 receptor
Symbol:	IL4R
UniProt ID:	P24394
Synonyms:	CD124; IL-4RA; IL4RA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)