

Product datasheet for **GA102411**

Human CD137 (TNFRSF9) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TNFRSF9 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 9 |
| Symbol: | TNFRSF9 |
| UniProt ID: | Q07011 |
| Synonyms: | 4-1BB; CD137; CDw137; ILA |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)