

## Product datasheet for **GA102445**

### Human EIF3E activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	EIF3E CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit E
Symbol:	EIF3E
UniProt ID:	<u><a href="#">P60228</a></u>
Synonyms:	eIF3-p46; EIF3-P48; EIF3S6; INT6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)