

Product datasheet for **GA102447**

Human IPP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IPP CRISPRa kit - CRISPR gene activation of human intracisternal A particle-promoted polypeptide
Symbol:	IPP
Components:	GA102447G1, IPP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102447G2, IPP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102447G3, IPP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y573</u>
Synonyms:	KLHL27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)