

Product datasheet for **GA102465**

Human CD41 (ITGA2B) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ITGA2B CRISPRa kit - CRISPR gene activation of human integrin subunit alpha 2b |
| Symbol: | ITGA2B |
| UniProt ID: | P08514 |
| Synonyms: | BDPLT2; BDPLT16; CD41; CD41B; GP2B; GPIIb; GT; GTA; HPA3; PPP1R93 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)