

Product datasheet for **GA102466**

Human Integrin alpha 3 (ITGA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITGA3 CRISPRa kit - CRISPR gene activation of human integrin subunit alpha 3
Symbol:	ITGA3
UniProt ID:	P26006
Synonyms:	CD49C; FRP-2; GAP-B3; GAPB3; ILNEB; MSK18; VCA-2; VL3A; VLA3a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)