

Product datasheet for **GA102468**

Human Integrin alpha 5 (ITGA5) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ITGA5 CRISPRa kit - CRISPR gene activation of human integrin subunit alpha 5 |
| Symbol: | ITGA5 |
| UniProt ID: | P08648 |
| Synonyms: | CD49e; FNRA; VLA-5; VLA5A |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)