

Product datasheet for **GA102488**

Human ITIH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITIH3 CRISPRa kit - CRISPR gene activation of human inter-alpha-trypsin inhibitor heavy chain 3
Symbol:	ITIH3
UniProt ID:	Q06033
Synonyms:	H3P; ITI-HC3; SHAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

