


Product datasheet for **GA102489**

Human ITIH4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITIH4 CRISPRa kit - CRISPR gene activation of human inter-alpha-trypsin inhibitor heavy chain 4
Symbol:	ITIH4
Components:	GA102489G1, ITIH4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102489G2, ITIH4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102489G3, ITIH4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14624</u>
Synonyms:	GP120; H4P; IHRP; ITI-HC4; ITIHL1; PK-120; PK120
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)