

Product datasheet for **GA102496**

Human IP3 receptor (ITPR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITPR1 CRISPRa kit - CRISPR gene activation of human inositol 1,4,5-trisphosphate receptor type 1
Symbol:	ITPR1
UniProt ID:	Q14643
Synonyms:	ACV; CLA4; INSP3R1; IP3R; IP3R1; PPP1R94; SCA15; SCA16; SCA29
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)