


## Product datasheet for **GA102498**

### Human ITPR3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ITPR3 CRISPRa kit - CRISPR gene activation of human inositol 1,4,5-trisphosphate receptor type 3
Symbol:	ITPR3
UniProt ID:	<a href="#">Q14573</a>
Synonyms:	IP3R; IP3R3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)