

Product datasheet for **GA102504**

Human JAK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	JAK3 CRISPRa kit - CRISPR gene activation of human Janus kinase 3
Symbol:	JAK3
UniProt ID:	<u>P52333</u>
Synonyms:	JAK-3; JAK3_HUMAN; JAKL; L-JAK; LJAK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)