

Product datasheet for **GA102513**

Human **KCNA1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 1
Symbol:	KCNA1
UniProt ID:	Q09470
Synonyms:	AEMK; EA1; HBK1; HUK1; KV1.1; MBK1; MK1; RBK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)