

Product datasheet for **GA102514**

Human Kv1.2 (KCNA2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 2
Symbol:	KCNA2
UniProt ID:	P16389
Synonyms:	EIEE32; HBK5; HK4; HUKIV; KV1.2; MK2; NGK1; RBK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)