

Product datasheet for **GA102515**

Human **KCNA3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 3
Symbol:	KCNA3
UniProt ID:	P22001
Synonyms:	HGK5; HLK3; HPCN3; HUKIII; KV1.3; MK3; PCN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)