

Product datasheet for **GA102516**

Human Kv1.4 (KCNA4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 4
Symbol:	KCNA4
UniProt ID:	P22459
Synonyms:	HBK4; HK1; HPCN2; HUKII; KCNA4L; KCNA8; KV1.4; PCN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)