

## Product datasheet for **GA102517**

### Human **KCNA5** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNA5 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 5
Symbol:	KCNA5
UniProt ID:	<a href="#">P22460</a>
Synonyms:	ATFB7; HCK1; HK2; HPCN1; KV1.5; PCN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)