

Product datasheet for **GA102518**

Human Kv1.6 (KCNA6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA6 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 6
Symbol:	KCNA6
UniProt ID:	P17658
Synonyms:	HBK2; KV1.6; PPP1R96
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)