


Product datasheet for **GA102519**

Human **KCNA7** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA7 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 7
Symbol:	KCNA7
UniProt ID:	Q96RP8
Synonyms:	HAK6; KV1.7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)