

Product datasheet for **GA102520**

Human Kv1.8 (KCNA10) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNA10 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member 10
Symbol:	KCNA10
UniProt ID:	Q16322
Synonyms:	Kcn1; Kv1.8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)