

Product datasheet for **GA102522**

Human **KCNC1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNC1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily C member 1
Symbol:	KCNC1
UniProt ID:	P48547
Synonyms:	EPM7; KV3.1; KV4; NGK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)