

Product datasheet for **GA102523**

Human Kv3.2 (KCNC2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNC2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily C member 2
Symbol:	KCNC2
UniProt ID:	<u>Q96PR1</u>
Synonyms:	KV3.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)