

Product datasheet for **GA102525**

Human KCNC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNC4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily C member 4
Symbol:	KCNC4
UniProt ID:	Q03721
Synonyms:	C1orf30; HKSHIIIC; KSHIIIC; KV3.4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)