

Product datasheet for **GA102526**

Human **KCND1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCND1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily D member 1
Symbol:	KCND1
UniProt ID:	<u>Q9NSA2</u>
Synonyms:	KV4.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)