


Product datasheet for **GA102528**

Human Kv4.3 (KCND3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCND3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily D member 3
Symbol:	KCND3
UniProt ID:	Q9UK17
Synonyms:	KCND3L; KCND3S; KSHIVB; KV4.3; SCA19; SCA22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)