

Product datasheet for **GA102529**

Human **KCNE1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNE1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily E regulatory subunit 1
Symbol:	KCNE1
UniProt ID:	P15382
Synonyms:	ISK; JLNS; JLNS2; LQT2/5; LQT5; MinK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)