

Product datasheet for **GA102530**

Human **KCNF1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNF1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily F member 1
Symbol:	KCNF1
UniProt ID:	<u>Q9H3M0</u>
Synonyms:	IK8; KCNF; KH1; KV5.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)