

Product datasheet for **GA102532**

Human **KCNH1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNH1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 1
Symbol:	KCNH1
UniProt ID:	O95259
Synonyms:	EAG; EAG1; h-eag; hEAG1; Kv10.1; TMBTS; ZLS1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)