

## Product datasheet for **GA102533**

### Human **KCNH2** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNH2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 2
Symbol:	KCNH2
UniProt ID:	<a href="#">Q12809</a>
Synonyms:	ERG-1; ERG1; H-ERG; HERG; HERG1; Kv11.1; LQT2; SQT1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)