

Product datasheet for **GA102536**

Human GIRK1 (KCNJ3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 3
Symbol:	KCNJ3
UniProt ID:	P48549
Synonyms:	GIRK1; KGA; KIR3.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)