

## Product datasheet for **GA102537**

### Human KIR2.3 (KCNJ4) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 4
Symbol:	KCNJ4
UniProt ID:	<u>P48050</u>
Synonyms:	HIR; HIRK2; HRK1; IRK-3; IRK3; Kir2.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)