

Product datasheet for **GA102538**

Human **KCNJ5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ5 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 5
Symbol:	KCNJ5
UniProt ID:	<u>P48544</u>
Synonyms:	CIR; GIRK4; KATP1; KIR3.4; LQT13
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)