

Product datasheet for **GA102539**

Human **GIRK2 (KCNJ6)** activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | KCNJ6 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 6 |
| Symbol: | KCNJ6 |
| UniProt ID: | P48051 |
| Synonyms: | BIR1; GIRK-2; GIRK2; hiGIRK2; KATP-2; KATP2; KCNJ7; KIR3.2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)