

Product datasheet for **GA102542**

Human Kir4.1 (KCNJ10) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ10 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 10
Symbol:	KCNJ10
UniProt ID:	P78508
Synonyms:	BIRK-10; KCNJ13-PEN; KIR1.2; KIR4.1; SESAME
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)