

Product datasheet for **GA102544**

Human Kir2.2 (KCNJ12) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ12 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 12
Symbol:	KCNJ12
UniProt ID:	Q14500
Synonyms:	hIRK; hIRK1; hkir2.2x; IRK-2; IRK2; kcnj12x; KCNJN1; Kir2.2; Kir2.2v
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)