

Product datasheet for **GA102545**

Human Kir7.1 (KCNJ13) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNJ13 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily J member 13
Symbol:	KCNJ13
UniProt ID:	<u>O60928</u>
Synonyms:	KIR1.4; KIR7.1; LCA16; SVD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)