

## Product datasheet for **GA102549**

### Human **KCNK1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNK1 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 1
Symbol:	KCNK1
UniProt ID:	<u><a href="#">O00180</a></u>
Synonyms:	DPK; HOHO; K2P1; K2p1.1; KCNO1; TWIK-1; TWIK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)