

Product datasheet for **GA102550**

Human KCNK2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | KCNK2 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 2 |
| Symbol: | KCNK2 |
| UniProt ID: | <u>O95069</u> |
| Synonyms: | hTREK-1c; hTREK-1e; K2p2.1; TPKC1; TREK; TREK-1; TREK1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)