

Product datasheet for **GA102551**

Human KCNK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK3 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 3
Symbol:	KCNK3
UniProt ID:	O14649
Synonyms:	K2p3.1; OAT1; PPH4; TASK; TASK-1; TBAK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)