

## Product datasheet for **GA102556**

### Human **KCNN3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNN3 CRISPRa kit - CRISPR gene activation of human potassium calcium-activated channel subfamily N member 3
Symbol:	KCNN3
UniProt ID:	<a href="#">Q9UGI6</a>
Synonyms:	hSK3; KCa2.3; SK3; SKCA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)