

## Product datasheet for **GA102558**

### Human **KCNQ1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNQ1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily Q member 1
Symbol:	KCNQ1
UniProt ID:	<a href="#">P51787</a>
Synonyms:	ATFB1; ATFB3; JLNS1; KCNA8; KCNA9; Kv1.9; Kv7.1; KVLQT1; LQT; LQT1; RWS; SQT2; WRS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)