

Product datasheet for **GA102559**

Human **KCNQ2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNQ2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily Q member 2
Symbol:	KCNQ2
UniProt ID:	<u>O43526</u>
Synonyms:	BFNC; BFNS1; EBN; EBN1; EIEE7; ENB1; HNSPC; KCNA11; KV7.2; KVEBN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)