

Product datasheet for **GA102560**

Human **KCNQ3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNQ3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily Q member 3
Symbol:	KCNQ3
UniProt ID:	O43525
Synonyms:	BFNC2; EBN2; KV7.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)