

Product datasheet for **GA102563**

Human **KCNS3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNS3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel modifier subfamily S member 3
Symbol:	KCNS3
UniProt ID:	Q9BQ31
Synonyms:	KV9.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)