

Product datasheet for **GA102572**

Human KIFC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIFC3 CRISPRa kit - CRISPR gene activation of human kinesin family member C3
Symbol:	KIFC3
Components:	GA102572G1, KIFC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102572G2, KIFC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102572G3, KIFC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BVG8</u>
Synonyms:	DKFZp686D23201; FLJ34694
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)