

Product datasheet for **GA102575**

Human KIR2DL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIR2DL3 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, two Ig domains and long cytoplasmic tail 3
Symbol:	KIR2DL3
UniProt ID:	P43628
Synonyms:	CD158b; CD158B2; GL183; KIR-023GB; KIR-K7b; KIR-K7c; KIRCL23; MGC129943; NKAT; NKAT2; NKAT2A; NKAT2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)