

Product datasheet for **GA102576**

Human KIR2DL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIR2DL4 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, two Ig domains and long cytoplasmic tail 4
Symbol:	KIR2DL4
UniProt ID:	<u>Q99706</u>
Synonyms:	CD158D; G9P; KIR-103AS; KIR103; KIR103AS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)